Organ Recovery
s y s t e m s

Chicago
2570 E. Devon Avenue
Des Plaines, IL 60018
tel 847.824.2600
tell free in US 866.682.4800
www.organ-recovery.com

Brussels
Da Vincilaan Box 6
1935 Zaventem
Belgium
tel 32 (0)2 715.0000
www.organ-recovery.com

European Authorized Representative
Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands

Representante Autorizado para Europa
Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands

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An ISO 13485 certified company
© Organ Recovery Systems, Inc. Printed in the USA

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Organ Recovery
s y s t e m s

LifePort™
Kidney Transporter Sterile Drape
LifePort™ disposable sterile drape for use only with LifePort™ perfusion equipment

Manufactured for:
Organ Recovery Systems
2570 E. Devon Avenue
Des Plaines, IL 60018
USA

847.824.2600

www.organ-recovery.com

Indications: Sterile Drape for LifePort™ Kidney Transporter is used to aid in maintaining aseptic conditions within the perfusion circuit. Follow hospital or physician-prescribed aseptic protocols.

Warnings:

 Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

 Before the operator begins to apply the sterile drape, the outer cassette lid should be removed. The outer surfaces of the outer cassette lid are not considered part of a sterile field; however, the inner surfaces should be handled as sterile.

 The application of the sterile drape is completed by a gowned and gloved operator. The inner surfaces of the LifePort™ Kidney Perfusion circuit are part of a sterile field. Use proper aseptic procedure in handling the drape and these components.

 The sterile drape is intended for a single use. After use, the drape should be disposed in accordance with local guidelines for biohazardous waste.

Instructions for use:

1. Carefully place the drape onto the kidney cassette, lining up the drape gasket and the well of the outer cassette lid. An orientation guide has been provided on the drape. The arrow on the guide should point towards the pump deck.
2. Unfold the drape, covering the transporter and exposing the inner cassette lid. The drape should unfold in the following order: right, left, front, back. The drape gasket should fit securely into the well of the outer cassette lid.
3. Unlatch and remove the inner cassette lid. Place the lid face-down onto a sterile surface.
4. Perform necessary procedure inside the perfusion circuit.
5. Replace inner cassette lid before removing drape.
6. Remove drape carefully to maintain the aseptic conditions of the inner cassette lid either cutting the drape and allowing it to fall away from the cassette or lifting the drape up and away from the field.
7. Replace outer cassette lid.

Storage condition: For best results, store at room temperature. Avoid excessive heat and humidity. Sterile unless package is damaged or open.

Contact: For technical assistance or to re-order please contact Organ Recovery Systems at 866.682.4800 (toll free) or 847.824.2600.

Caution: Federal and international law restricts this device to sale by or on the order of a physician or licensed practitioner.

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