Organ Recovery
systems

Chicago
2570 E. Devon Avenue
Des Plaines, IL 60018
tel 847.824.2600
toll free in US 866.682.4800
www.organ-recovery.com

Brussels
Da Vincilaan Box 6
1935 Zaventem
Belgium
tel 32 (0)2 715.0000
www.organ-recovery.com

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
DEUTCH

LifePort® Kidney Transporter Disposable Cannula

SealRing™

CAN-0720
CAN-1035

Indication: The SealRing™ cannula is used when the vessel to be perfused terminates with an aortic patch or similar.

1. Use standard aseptic procedures to introduce the cannula onto the sterile field.
2. Slide the aortic patch through the center of the cannula ring and lay patch flat, such that the tissue covers the entire SealRing. If necessary, instruments may be used to temporarily hold the tissue in place until the cannula is secure.
3. Lower the upper portion of the cannula securing the tissue between the two halves.
4. Wrap each strap around both cannula halves and secure the strap to the post.
5. Snap the cannula into the cannula mount.
6. Adjust the height and rotation of the cannula mount to position the vessel. Visually inspect the vessel, ensuring there are no twists or occlusions.
7. Connect the infuse line and remove the end cap.
8. Follow LifePort® perfusion equipment instructions to remove air from the infuse line and the cannulated vessel.
9. Replace end cap.
10. Begin infusion in accordance with perfusion equipment instructions and surgeon-specified protocol.
11. Visually check for leaks. Leaks may be sealed by repositioning the tissue in the cannula or by adjusting the tension on the cannula straps.

Coupler


CAN-1000

Indication: The Coupler is used to connect two cannula when multiple vessels must be perfused. Indication: The Capler is used to connect two cannula when multiple vessels must be perfused.

1. Use standard aseptic procedures to introduce the coupler and required cannulas onto the sterile field.
2. Apply a SealRing™ or straight cannula to each of the vessels to be perfused.
3. Identify a main vessel.
4. Place the organ in the kidney cradle and snap the cannula connected to the main vessel into the cannula mount.
5. Replace the end cap with one end of the coupler.
6. Attach the second end of the coupler to the infuse port of the next cannula.
7. Repeat 5 and 6 for each cannula being used.
8. Connect the infuse line to the main vessel cannula and remove the end cap of the last cannula in the chain.
9. Verify the coupler tubing or vessels are neither kinked nor occluded.
10. Follow LifePort® perfusion equipment instructions to remove air throughout the infuse line, cannulas, coupler tubing, and vessels.
11. Replace end cap.
12. Begin infusion in accordance with perfusion equipment instructions and surgeon-specified protocol.
13. Visually check for and resolve all leaks.

Straight

CAN-0003
CAN-0005
CAN-0008

Indication: the straight cannula is used when the vessel to be perfused terminates without a patch or when intimal lining damage is not a concern.

1. Use standard aseptic procedures to introduce the cannula onto the sterile drape.
2. Place the cannula tip into the vessel.
3. Secure the vessel in place using silk ties, vessel loops,

or another appropriate material. A groove in the tip is provided for positioning.

4. Snap the cannula into the cannula mount.
5. Adjust the height and rotation of the cannula mount to position the vessel. Visually inspect the vessel, ensuring there are no twists or occlusions.
6. Connect the infuse line and remove the end cap.
7. Follow LifePort® perfusion equipment instructions to remove air from the infuse line and the cannulated vessel.
8. Replace end cap.
9. Begin infusion in accordance with perfusion equipment instructions and surgeon-specified protocol.
10. Visually check for leaks. Leaks may be sealed by repositioning the tissue or increasing tension on ties or loops.

Storage Conditions: For best results, store at room temperature. Avoid excessive heat and humidity. Sterile unless package is damaged or open.

Contact: For technical assistance or to re-order please contact Organ Recovery Systems at 866.682.4800 (toll free) or 847.824.2600.

Caution: Federal and international law restricts this device to sale by or on the order of a physician or licensed practitioner.

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© 2006 Organ Recovery Systems, Inc. Printed in the USA